

Wane CORE_{XP}



Decoding the future



Conical Connection :

Ensures a secure fit with abutment, eliminates micro movements reducing bone loss at crestal level.




Reverse Conical Neck :

Preserves cortical bone, promotes alveolar bone remodeling.



Conical Shape with deep, sharp threads :

Gradual bone condensing effect, self drilling, self tapping, controlled and easy insertion in low density bone.



Apical Taper :

Quick insertion into the bone. Optimal anchorage on the apical portion.

Connection :

Internal hex conical connection with screw retained abutments.

Available in two platform switching of 2.8 & 3.4 mm for 3.25/3.5/3.9 & 4.3/5.0/5.5 platform respectively.

Diameters	3.25	3.5	3.9	4.3	5.0	5.5
------------------	------	-----	-----	-----	-----	-----

Lenghts	8.0	10.0	11.5	13.0	15.0
----------------	-----	------	------	------	------

Dimensions in mm

