



Decoding the future

Wane CORE





Conical Connection:

Ensures a Secure fit with abutment, eliminates micro movements reducing bone loss at crestal level.

Coronal Threads:

Enhances primary stability and preserves the crestal bone.

Body Threads:

In addition to SLA surface, dual thread design increases surface area which leads to higher bone to implant contact (BIC).

Apical Taper:

Enable to easily place the implant into ossteotomy thereby completely securing placement.

Connection:

Internal hex conical connection with screw retained abutments.

Available in two platform switching of 2.8 & 3.4 mm for 3.5/3.9 & 4.3/5.0 platform respectively.

Tanalita a a da a da a da a da a	Diameters	3.5	3.9	4.3	5.0	
Lenghts 8.0 10.0 11.5 13.0 15	Lenghts	8.0	10.0	11.5	13.0	15.0

Dimensions in mm

Prosthetics

Type of Abutment		Connection	Emergence	Gingival Height	Abt. Height
Hollow		3.25/3.5/3.9	3.5/4.0/5.0	0.75/1.5/2.5/3.5	5.5/7.0
	Abutment	4.3/5.0/5.5	4.0/5.0/6.0	0.75/1.5/2.5/3.5	5.5/7.0
Angled Abutment	3.25/3.5/3.9	3.5/4.0/5.0	0.75/1.5/2.5/3.5	6.4	
A	(15°/25°)		4.0/5.0/6.0	0.75/1.5/2.5/3.5	6.4
Ball Abutment	3.25/3.5/3.9	3.5	0.0/0.75/1.5/2.5/3.5	3.0	
	Ball Abutment	4.3/5.0/5.5	4.0	0.0/0.75/1.5/2.5/3.5	3.0
Ti Paga Ala	Ti-Base Abutment	3.25/3.5/3.9	Hex/Non Hex	0.75/1.5/2.5/3.5	-
71	II-Dase Aputilient	4.3/5.0/5.5	Hex/Non Hex	0.75/1.5/2.5/3.5	-
Castable	3.25/3.5/3.9	Hex/Non Hex	-	-	
VU.	Abutment	4.3/5.0/5.5	Hex/Non Hex	-	-
Straight Multi	3.25/3.5/3.9	4.8	0.75/1.5/2.5/3.5	2.2	
41	Abutment	4.3/5.0/5.5	4.8	0.75/1.5/2.5/3.5	2.2
Angled Multi Abutment	3.25/3.5/3.9	4.8	0.75/1.5/2.5/3.5	2.2	
HY	(17°/30°)	4.3/5.0/5.5	4.8	0.75/1.5/2.5/3.5	2.2
Gingival Former		3.25/3.5/3.9	3.5/4.0/5.0	0.75/1.5/2.5/3.5	3.0
	dingival rolline	4.3/5.0/5.5	4.0/5.0/6.0	0.75/1.5/2.5/3.5	3.0
Transfer Coping Open Tray : Available in : 3.25/3.5/3.9 4.3/5.0/5.5		Transfer Coping Open Tray Multi Abutment : Available in: 3.25/3.5/3.9 4.3/5.0/5.5		Lab Analog - Implant : Available in : 3.25/3.5/3.9 4.3/5.0/5.5	
Transfer Coping Close Tray : Available in : 3.25/3.5/3.9 4.3/5.0/5.5		Transfer Coping Close Tray Multi Abutment : Available in: 3.25/3.5/3.9 4.3/5.0/5.5		: Available in : 3.25/3.5/3.9 4.3/5.0/5.5	

Prosthetics: Wide Range of abutments for fixed, fixed-detachable and screw retain prosthetics.










Surgical Motor



Specifications

 Overload Protection 	• 30-50000 Rpm Speed
• 20:1 Reduction	• 5-55 Ncm Torque
• 10 Programmable memories	• 8 Gear Ratios
Dual Operation	 4 Water Speeds





Mobile / WhatsApp : +91 9321865230 / +91 7567379882

Website: www.veritasbio.com

Email: mail@veritasbio.com / export@veritasbio.com

